Taller 2 computacional

Lizeth Gomez

2026-04-12

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Se utilizó la herramienta **PrimerBlast**:

Moniliophthora roreri

The image shows the PrimerBlast web interface, which is used for designing primers. It is divided into several sections:

- PCR Template:** This section allows users to enter an accession, GI, or FASTA sequence. It also includes a "Range" selector and fields for "Forward primer" and "Reverse primer". A file upload option is available, with a file named "M_roreri.fasta" selected.
- Primer Parameters:** This section contains fields for "Use my own forward primer (5'->3' on plus strand)" and "Use my own reverse primer (5'->3' on minus strand)". It also includes a "PCR product size" range (Min: 100, Max: 500), a "# of primers to return" field (set to 10), and "Primer melting temperatures (T_m)" (Min: 57.0, Opt: 60.0, Max: 63.0). A "Primer Location Preference" checkbox is also present.
- Exon/intron selection:** This section includes a "Refseq mRNA sequence as PCR template input is required for options in the section" note. It has fields for "Exon junction span" (set to "No preference"), "Exon junction match" (Min 5' match: 7, Min 3' match: 4, Max 3' match: 8), and "Intron inclusion" (checkbox). A "Minimal and maximal number of bases that must anneal to exons at the 5' or 3' side of the junction" field is also present.
- Primer Pair Specificity Checking Parameters:** This section includes a "Note: Parameter values that differ from the default are highlighted in yellow". It has a "Specificity check" section with a "Search mode" dropdown (set to "Automatic") and a "Database" dropdown (set to "Refseq representative genomes"). It also includes a "Max target amplicon size" field (set to 4000) and a "Allow splice variants" checkbox.

Moniliophthora perniciosa

Retrieve recent results Publication Tips for finding specific primers

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [?](#) [Clear](#)

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Primer Location Preference ☐ Prefer primers at 3' side of the template [?](#)

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section [?](#)

Exon junction span No preference [?](#)

Exon junction match Min 5' match Min 3' match Max 3' match

Minimal and maximal number of bases that must anneal to exons at the 5' or 3' side of the junction [?](#)

Primer Pair Specificity Checking Parameters

Note: Parameter values that differ from the default are highlighted in yellow

Specificity check ☒ Enable search for primer pairs specific to the intended PCR template [?](#)

Search mode Automatic [?](#)

Database Refseq representative genomes [?](#)

Exclusion ☐ Exclude predicted Refseq transcripts (accession with XM, XR prefix) ☐ Exclude uncultured/environmental sample sequences [?](#)

Organism Marasmiaceae (taxid:654128) [Add organism](#)

Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type. [?](#)

Entrez query (optional) [?](#)

Primer specificity stringency Primer must have at least 2 total mismatches to unintended targets, including at least 2 mismatches within the last 5 bps at the 3' end. [?](#)

Ignore targets that have 6 or more mismatches to the primer. [?](#)

Max target amplicon size 4000 [?](#)

Allow splice variants ☐ Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input) [?](#)

[Get Primers](#) ☒ Show results in a new window ☒ Use new graphic view [?](#)

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Primer Location Preference ☐ Prefer primers at 3' side of the template [?](#)

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section [?](#)

Exon junction span [?](#)

Exon junction match

Minimal and maximal number of bases that must anneal to exons at the 5' or 3' side of the junction [?](#)

Intron inclusion ☐ Primer pair must be separated by at least one intron on the corresponding genomic DNA [?](#)

Intron length range [?](#)

Primer Pair Specificity Checking Parameters

Note: Parameter values that differ from the default are highlighted in yellow

Specificity check ☒ Enable search for primer pairs specific to the intended PCR template [?](#)

Search mode [?](#)

Database [?](#)

Exclude predicted Refseq transcripts (accession with XM, XR prefix) ☐ Exclude uncultured/environmental sample sequences [?](#)

Organism [Add organism](#)

Primer-BLAST JOB ID:amC3cU_3946ds14s-3q2K7-yANPTng0

Primer-BLAST Results [?](#)

Input PCR template: none [Edit search](#) [New](#)

Specificity of primers: Target templates were found in selected database: Core nucleotide BLAST database (Organism limited to Fungi) [? Search history](#)

Other reports [?](#)

Detailed primer reports

[Download primer pairs](#)

Primer pair 1

Sequence (5'→3')	Length	Tm	GC%	Self complementarity	Self 3' complementarity
Forward primer	20	58.71	50.00	4.00	2.00
Reverse primer	20	60.18	50.00	2.00	2.00

Products on target templates

>CC2919137.1 Monilophthora perniciosa isolate CPB17 18S ribosomal RNA gene, internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene, complete sequence; and internal transcribed spacer 2, partial sequence

product length = 413

Forward primer	Reverse primer	Template
1	1	1
15	15	24
Reverse primer	1	28
Template	429	486

M. perniciosa

Primer Parameters

Use my own forward primer (5'→3' on plus strand) [?](#) [Clear](#)

Use my own reverse primer (5'→3' on minus strand) [?](#) [Clear](#)

PCR product size

of primers to return

Primer melting temperatures (T_m) [?](#)

Primer Location Preference ☐ Prefer primers at 3' side of the template [?](#)

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section [?](#)

Exon junction span [?](#)

Exon junction match

Minimal and maximal number of bases that must anneal to exons at the 5' or 3' side of the junction [?](#)

Intron inclusion ☐ Primer pair must be separated by at least one intron on the corresponding genomic DNA [?](#)

Intron length range [?](#)

Primer Pair Specificity Checking Parameters

Note: Parameter values that differ from the default are highlighted in yellow

Specificity check ☒ Enable search for primer pairs specific to the intended PCR template [?](#)

Search mode [?](#)

Database [?](#)

Exclude predicted Refseq transcripts (accession with XM, XR prefix) ☐ Exclude uncultured/environmental sample sequences [?](#)

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